



User Guide

шшш.globescientific.com

## **TABLE OF CONTENTS**

Section	Description	Page #
1	Introduction	1
2	Package inclusion	1
3	Product Description	2
3.1	Adjustable Volume Pipette Range	2
3.2	Fixed volume pipette range	3
3.3	Multi channel pipette range	3
3.4	Setting the Volume	4
3.5	Volume indicator display example	4
4	Pipette operation	5
4.1	Good Pipetting Guidelines	5
4.2	Setting the volume	5
4.3	Loading tips	6
4.4	Optimum immersion depths	6
4.5	Forward pipetting	6
4.6	Reverse pipetting	7
4.7	Aspiration of sample	7
4.8	Dispensing sample	7
4.9	Ejection of tips	8
5	Calibration and adjustment	9
5.1	Calibration Requirements and Test Conditions	8
5.2	Calibration adjustment	8
5.3	Procedure to check calibration	9
5.3.1	Conversion of weight readings to volume	9
5.3.2	Calculation for in-accuracy (systematic error)	9
5.3.3	Calculation for imprecision (random error)	9
6	Maintenance and servicing	9
6.1	Disassembly	10
6.2	Disassembling the lower part	10
6.2.1	Disassembling the single channel pipette up to 1ml	10
6.2.2	Disassembling the single channel pipette: 5-10ml	11
6.3	Assembling the pipette	11
6.3.1	Assembling the single channel pipette up to 1ml	11
6.3.2	Assembling the single channel pipette: 5-10ml	11
6.4	Disassembling the multi channel pipette	11
6.4.1	Removing the lower assembly	11
6.4.2	Opening the lower assembly	12
6.4.3	Removing the channel	12
6.4.4	Fitting the channel	12
6.4.5	Assembling the lower assembly	13
6.4.6	Checking the function	12
6.5	List of spare parts - Single Channel Pipette	13
6.6	List of spare parts - Multichannel Pipette	18
6.7	Autoclaving	19
7	Trouble shooting guide	19

## 1. INTRODUCTION


Thank you for choosing the DiamondAPEX series of pipettors! Numerous advances in development have been implemented, including:

- Magnet assisted piston to aid in precise results
- Innovative spring and seal design for extreme reduction in plunger force required
- Corrosion resistant tip ejector mechanism with unique shockabsorbing design
- Volume lock
- Fully autoclavable
- Ergonomic design
- Easy in-house calibration
- Highly durable universal tip cone


## 2. PACKAGE CONTENTS

Des	cription	Quantity
•	DiamondAPEX Pipette	1
•	Certificate of conformity including calibration certifica	te 1
•	Warranty card	1
•	Product manual	1
•	Shelf-mounting stand	1
•	Calibration tool	1
•	Silicone grease	1
•	Sample pack of tips	1

## SINGLE CHANNEL ADJUSTABLE VOLUME



## **MULTICHANNEL ADJUSTABLE VOLUME**



## 3. PRODUCT DESCRIPTION

## 3.1 ADJUSTABLE VOLUME PIPETTE RANGE

CAT No.	Color Code	Volume Range [ μl ]	Incre- ments [ μl ]	Test Volume [ μl ]	Inac- curacy (±) %	Imprecision (±) %
				0.25	12	6
3351-1		0.1-2.5 μΙ	0.002	1.25	2.5	1.5
				2.5	2.5	0.7
				1	2.5	1.5
3351-10		0.5-10 μΙ	0.02	5	1.5	0.8
				10	1	0.4
				2	3	1.5
3351-20		2-20 µl	0.02	10	1.2	0.6
				20	0.9	0.3
				5	2	2
3351-50		5-50 μl	0.1	25	0.8	0.4
				50	0.6	0.3

			10	3	1
3351-100	10-100 μΙ	0.1	50	1	0.3
			100	0.8	0.2
			20	2.5	0.7
3351-200	20-200 μΙ	0.2	100	0.7	0.3
			200	0.6	0.2
			100	3	0.6
3351-1000	100-1000 μΙ	1	500	1	0.2
			1000	0.6	0.2
			500	2.4	0.6
3351-5000	500-5000 μΙ	10	2500	1.2	0.25
			5000	0.6	0.2
			1000	3	0.6
3351-10000	1000-10000 μΙ	20	5000	0.8	0.2
			10000	0.6	0.15

## 3.2 FIXED VOLUME PIPETTE RANGE

CAT No.	Color Code	Volume Range [ μl ]	Test Volume	Inaccuracy (±) %	Imprecision (±) %
	555		[µl]	, , ,	(-)
3352-2		2.5 μΙ	2.5	2	1.6
3352-5		5 μΙ	5	1.3	1.2
3352-10		10 μΙ	10	1.2	0.6
3352-20		20 μΙ	20	1	0.3
3352-25		25 μΙ	25	1	0.3
3352-50		50 μΙ	50	0.7	0.3
3352-100		100 μΙ	100	0.6	0.2
3352-200		200 μΙ	200	0.6	0.2
3352-250		250 μΙ	250	0.6	0.3
3352-500		500 μΙ	500	0.6	0.2
3352-1000		1000 μΙ	1000	0.6	0.2
3352-2000		2000 μΙ	2000	0.3	0.15
3352-5000		5000 μΙ	5000	0.3	0.15
3352-10000		10000 μΙ	10000	0.6	0.2

# 3.3 MULTI-CHANNEL PIPETTE RANGE 8 Channels & 12 Channels

CAT No.	Color Code	Volume Range [ μl ]	Incre- ments [ μl ]	Test Volume [ μl ]	Inac- curacy (±) %	Imprecision (±) %
				1	8	5
3354-10 3355-10		0.5-10 µl	0.02	5	4	2
				10	2	1

			2	7	3
3354-20 3355-20	2-20 μl	0.02	10	3	2
			20	2	1.6
			5	3	2
3354-50 3355-50	5-50 μl	0.1	25	1.5	1
			50	1	0.7
			10	3	2
3354-100 3355-100	10-100 μΙ	0.1	50	1	0.8
3333 230			100	0.8	0.3
			20	5	1.4
3354-200 3355-200	20-200 μΙ	0.2	100	1	0.4
			200	0.7	0.25
			30	3	1
3354-300 3355-300	30-300 μΙ	0.2	150	1	0.5
3333 300			300	0.6	0.3

## 3.4 SETTING THE VOLUME


Delivery volume is clearly indicated in the volume display seen in the main body of the pipettor.

In variable volume models, the bottom volume wheel includes a small increment scale for precise setpoint and delivery capabilities.

Some variable volume pipettors include one or two decimal places in the setpoint volume wheels. This is indicated by the use of a black horizontal line as seen in the following examples section.



## 3.5 VOLUME INDICATOR DISPLAY EXAMPLES



## 4 PIPETTE OPERATION

## 4.1 GOOD PIPETTING GUIDELINES

- Check the tip cone to make sure it is clean
- Make sure the pipette tip is firmly attached to the tip cone
- Ensure that the temperature of the pipette tip, pipettor and sample are at equilibrium
- Always press and release the plunger slowly and smoothly
- Pre-rinse the pipette tip before aspirating the sample by filling and emptying the pipette tip five (5) times. This is important when dispensing samples that have a viscosity and density different from water and for volatile solvents
- Do not pre-rinse the tip when pipetting samples with temperatures different from the current ambient temperature. Be sure to change the pipette tip after each pipetting.
- When aspirating, hold the pipettor in an upright position and keep the pipette tip at a constant depth below the surface of the sample liquid.
- After pipetting acids or other corrosive liquids, remove the tip cone and and rinse the piston, o-ring and seal with distilled water
- Do not pipette samples at a temperature above 70°C or below 4°C

Make sure that the sample never enters the tip cone. To prevent this:

- Never lay the pipettor on it's side when there is liquid in the pipette tip
- Always press and release the plunger slowly and smoothly
- Never turn the pipettor upside-down

#### 4.2 SETTING THE VOLUME

 To set the volume, turn the volume lock lever to the 'unlock' position so that the volume can be adjusted to the desired setpoint within the permitted volume range









- To decrease the volume, turn the plunger knob clockwise
- To increase the volume, turn the plunger knob counterclockwise
- Turn the volume lock lever to the 'lock' position to avoid accidentally changing the volume during pipetting



The locking mechanism ensures that the volume adjustment knob remains at the desired setpoint while aspirating or dispensing sample liquids. Any effort to rotate the volume adjustment knob with the locking mechanism engaged will damage the locking mechanism and void the warranty.

A Setting the volume beyond the allowable volume range is not permitted. Using excessive force to turn the volume adjustment knob outside the permitted range will jam the mechanism, damage the pipette and void the warranty.

#### 4.3 LOADING TIPS


Globe Scientific recommends using authentic Globe Scientific pipette tips to ensure a perfect fit for maximum precision.

- Use the correct volume pipette tip according to the volume range of the pipettor.
- Be sure to check the tip cone to ensure it is clean.
- Firmly press the pipette tip on to the tip cone to ensure an airtight seal.

## 4.4 OPTIMUM IMMERSION DEPTHS

Volume	Immersion depth
0.1 μl - 1 μl	1 mm
1 μl - 100 μl	2 - 3 mm
100 μΙ - 1000 μΙ	2 - 4 mm
1ml - 10ml	3 - 5 mm

Tip immersion depth is critical and should not be exceeded or the volume may become inaccurate. Tip angle is also important; the pipettor should be used only at an angle between zero and twenty degrees (0°-20°):



## 4.5 FORWARD PIPETTING

- 1. Press the plunger until the first stop.
- Dip the pipette tip under the surface of the sample liquid in the reagent reservoir and then slowly release the plunger. Withdraw the pipette tip from the liquid, touching it against the edge of the reagent reservoir to remove any excess liquid or droplets.

- Dispense the sample liquid by gently pressing the plunger to the first stop. After a delay of approximately one second, continue to press the plunger fully to the second stop, effectively emptying the pipette tip.
- Release the plunger fully (to the 'ready' position). Change the pipette tip and continue pipetting as necessary.



#### 4.6 REVERSE PIPETTING

The reverse pipetting technique is suitable for dispensing liquids that have high viscosity or a tendency to foam easily. This technique is also recommended for dispensing very small volumes.

- 1. Press the plunger fully to the second stop.
- Dip the pipette tip under the surface of the sample liquid in the reagent reservoir and then slowly release the plunger. This step will fill the tip. Withdraw the pipette tip from the sample liquid and touch it against the edge of the reagent reservoir to remove any excess sample liquid.
- 3. Dispense the preset volume of sample liquid by gently pressing the plunger to the first stop. Hold the plunger at the first stop. Some of the sample liquid Second Stop will remain in the pipette tip and it should not be dispensed.
- The remaining sample liquid should either be discarded with the pipette tip or should be dispensed back into the reagent reservoir.



## 4.7 ASPIRATION OF SAMPLE

- Hold the pipettor vertically and press the plunger to its first stop. Place the pipette tip into the sample liquid at proper angle and depth (please refer to OPTIMUM IMMERSION DEPTHS) then gently release your thumb pressure on the plunger. The plunger will then move upward, aspirating the sample.
- Pause for about 1 second (longer for macro-volume pipettors) to ensure that the full volume of sample is drawn into the pipette tip.
- Withdraw the pipette tip from the sample liquid. If any liquid remains on the outer surface of the pipette tip, touch it off carefully in a lint-free tissue, taking care not to touch the pipette tip orifice.

#### 4.8 DISPENSING SAMPLE

- 1. Place the pipette tip against the reagent reservoir wall to avoid any bubbles or splashing of the sample liquid out of the reagent reservoir.
- Press the plunger slowly past the first stop and second stop for complete blow out of the sample liquid. For viscous samples, it is recommended to wait for few seconds until the sample fluid passes out to the reservoir.
- 3. Pull the pipette tip along the wall of the reagent reservoir and

- release the plunger slowly.
- 4. Discard the pipette tip to avoid any sample carryover / cross-contamination. It is highly recommended to replace the pipette tip prior to repeating the pipetting cycle.

## 4.9 EJECTION TIPS

The pipette tip ejector should be pressed downwards firmly and fully to ensure proper tip ejection. Pipette tips should be ejected and disposed of in a suitable waste container.

## 5 CALIBRATION AND ADJUSTMENT

All pipettes have been quality tested according to ISO8655-6. The quality control process according to ISO 8655-6 involves gravimetric testing of each pipette with double distilled water.

All pipettes are calibrated in an ISO/IEC 17025 accredited laboratory. Each Pipette is calibrated, inspected and validated by qualified technicians according to defined quality systems.

## 5.1 CALIBRATION REQUIREMENTS AND TEST CONDITIONS

An analytical balance must be used. The balance selection depends upon the selected model of the pipettor and the sensitivity of the balance reading.

Test liquid: Water, distilled or deionized, grade 3 water conforming to ISO3696. Calibration should be performed in a draft-free room with the temperature of the water, pipettor and air between  $15^{\circ}$ C to  $30^{\circ}$ C (stable within  $\pm 0.5^{\circ}$ C).

The relative humidity must be above 50%, especially with calibration volumes under 50µl to help reduce evaporation loss.

Special accessories for the analytical balance such as an evaporation trap are recommended for calibration volumes under  $50\mu$ l.

## 5.2 CALIBRATION AND ADJUSTMENT

- 1. Calibration adjustment is performed with the calibration tool provided with the pipettor.
- 2. Rotate the volume locking lever to the "lock" position so the volume setting mechanism is locked and able to turn the calibration screw.





3. Remove the volume adjustment knob by pulling it upwards.





4. Place the calibration tool into the calibration grooves.



- 5. Turn the calibration tool counterclockwise to increase volume and clockwise to decrease volume.
- 6. After adjustment, check the calibration according to PROCEDURE TO CHECK CALIBRATION in the next section.

⚠ Depending upon use, we recommend a calibration check frequency of every six months. This schedule is only a general recommendation and should be adjusted according to your institution's use and requirements.

## 5.3 PROCEDURE TO CHECK CALIBRATION

The pipettor is checked at maximum volume, at 50 % of maximum volume and at minimum or 10% of maximum volume, whichever is higher.

- A new tip is first pre-wetted 3-5 times and a series of ten pipetting is done at each volume.
- Use of forward pipetting technique is recommended.
- Calculate the inaccuracy and imprecision for all three volumes as per EN ISO 8655-6 standards on the basis of the following calculation

#### 5.3.1 CONVERSION OF WEIGHT READINGS TO VOLUME

Mean Volume  $\overline{V} = \overline{X} \cdot Z$ Mean Weight  $\overline{X} = \underbrace{\sum Xi}_{n}$ 

Xi = Balance Reading
n = number of reading
Z = Conversion factor

[example Z=1.0040 μl/mg at 25°c and 1013 hPa]

## 5.3.2 CALCULATION FOR IN-ACCURACY (SYSTEMATIC ERROR)

 $A\% = \frac{\overline{V} - V_0}{V_0} \cdot 100$ 

 $\overline{V}$  = Mean Value

Vo = Particular volume at which readings are taken

## 5.3.3 CALCULATION FOR IMPRECISION (RANDOM ERROR)

$$S = \sqrt{\frac{\sum_{i=1}^{n} (Vi - \overline{V})^2}{n - 1}}$$

$$CV\% = \frac{100 \cdot S}{\overline{V}}$$

S = Standard Deviation

V = Mean Value

n = number of readings

Compare the results to the limits in the earlier tables (Page# 4, 5)

## 6 MAINTENANCE AND SERVICE

When the pipettor is not in use it should be stored in an upright position. The pipettor should be inspected prior to use each day for

any dust or contamination on outside surfaces. Special attention should be given to the tip cone. No solvent other than isopropanol should be used to clean the pipette. If the pipettor is used daily, an internal parts inspection should be performed every three (3) months.

## 6.1 DISASSEMBLY

The servicing procedure starts with the disassembly of the pipettor. Please refer to the spare parts lists for better understanding of the components.


#### 6.2 DISASSEMBLING THE LOWER PART

- 1. Press tip ejector pusher completely down and hold.
- 2. Pull down the tip ejector and release the tip ejector pusher.



## 6.2.1 DISASSEMBLING THE SINGLE-CHANNEL PIPETTE UP TO 1ml

- 1. Unscrew the coupler and remove tip cone.
- 2. Remove the lower part and pull out the piston and other parts from tip cone.
- Remember to keep all parts in order for reassembly.
- Clean the piston, the piston spring, seal and the o-ring with isopropanol and lint free tissue.
   Allow them to dry.
- Check the tip cone for foreign articles and remove, if any. Grease the cleaned parts with the approved lubricant provided with each pipettor.



### 6.2.2 DISASSEMBLING THE SINGLE-CHANNEL PIPETTE: 5 - 10ml

- 1. As per earlier process remove the ejector by pulling it after pressing the tip ejector pusher completely down.
- Unscrew the coupler from tip cone and pull out the piston and other parts from tip cone.
- Remember to keep all parts in order for reassembly.
- Clean the piston, the piston spring, seal and the o-ring with isopropanol and lint free tissue. Allow them to dry.
- Check the tip cone for foreign particles and remove, if any. Grease the cleaned parts with the approved lubricant provided with each pipette.



## 6.3 ASSEMBLING THE PIPETTE

## 6.3.1 ASSEMBLING THE SINGLE-CHANNEL PIPETTE UP TO 1ml

- 1. Carefully insert the piston into the tip cone.
- 2. Press on piston from above and check for free movement. The piston must be able to move freely without resistance.
- Reconnect the tip cone to the main body by screwing into the threaded section.
- 4. Reinstall the tip ejection collar.

## 6.3.2 ASSEMBLING THE SINGLE-CHANNEL PIPETTE: 5 - 10ml

- 1. Keep retainer on piston and spring. Press the spring to fit with piston.
- 2. Carefully insert the piston into the tip cone.
- 3. Press on piston from above and check for free movement. The piston must be able to move freely without resistance.
- Reconnect the tip cone to the main body by screwing into the threaded section.
- 5. Reinstall the tip ejection collar.

# 6.4 DISASSEMBLING THE MULTI-CHANNEL PIPETTE 6.4.1 REMOVING THE LOWER ASSEMBLY

- 1. Press the tip ejector pusher completely and hold it while unscrewing the coupler from upper part of pipette.
- 2. Remove the lower assembly.



## 6.4.2 OPENING THE LOWER ASSEMBLY

- 1. Unscrew the two small screws from back side and safely keep them.(Refer to image -1)
- 2. Press and push up from side to open the front cover. (Refer to image -2)







#### 6.4.3 REMOVING THE CHANNEL

- 1. Slightly push the spring and pull up tip cone to remove it from the lower rail. (image -3)
- 2. Carefully release the piston from the upper rail and move upward to remove it. (image -4)





## 6.4.4 FITTING THE CHANNEL

- 1. Insert the spring with cylinder into the center rail.
- 2. Insert the piston into the cylinder and fit into upper rail.
- 3. Compress the spring with the cylinder and insert the cylinder into the lower rail.

## 6.4.5 ASSEMBLING THE LOWER ASSEMBLY

- 1. Attach the front cover and screw it in.
- 2. Press the pusher completely and hold it while screwing the coupler nut to body.

#### 6.4.6 CHECKING THE FUNCTION

It is to ensure that the pipettor has been assembled correctly.

- Carry out a gravimetric test for the systematic and random error.

## 6.5 LIST OF SPARE PARTS - SINGLE CHANNEL PIPETTE

Group 2

Group 1

Adjustable volume 0.1 -Adjustable volume 0.5-10μl 2.5µl Fix volume 2.5µl Fix volume 10µl DP 030 DP 031 DP 0269 DP 0269 DP 036-DP 037 DP 0261 DP 0262 DP 0159 DP 0159 SP 0158 SP 0160-DP 022-DP 022


Group 3 Adjustable volume 2-20µl

Fix volume 20μΙ Group 4

Adjustable volume 5-50µl


Fix volume

25μΙ 50µl



**Group 5** Adjustable volume 10-100μl Fix volume 100μl

Group 6
Adjustable volume 20200μl
Fix volume 200μl



## Group 7

Adjustable volume  $100 - 1000 \mu l$ 

Fix volume 250µl

500μΙ


1000μΙ



**Group 8** Group 9 Adjustable volume 0.5-5ml Adjustable volume 1-10ml Fix volume 2ml Fix volume 10ml 5ml DP 021 DP 020-DP 0268 DP 0268 DP 0273 DP 0272 SP 0168 ----SP 0168 **←**SP 0044 SP 0060 -SP 0162 RP 0156 SP 0108 SP 0110 -



## 6.6 LIST OF SPARE PARTS - MULTICHANNEL PIPETTE



Part	Cat#							
	3354-10	3354-20	3354-50	3354-100	3354-200	3354-300		
	3355-10	3355-20	3355-50	3355-100	3355-200	3355-300		
	(0.5-10µl)	(2-20µl)	(5-50µl)	(10-100µI)	(20-200µl)	(30-300µl)		
А	DPM 165	DPM 165	DPM 165	DPM 165	DPM 165	DPM 165		
В	DPM 134	DPM 135	DPM 136	DPM 137	DPM 138	DPM 139		
С	DPM 123	DPM 123	DPM 124	DPM 125	DPM 126	DPM 127		
D	DPM 276	DPM 276	DPM 277	DPM 277	DPM 277	DPM 277		
E	DPM 364	DPM 364	DPM 365	DPM 366	DPM 367	DPM 368		
F	SPM 152	SPM 154	SPM 101	SPM 155	SPM 150	SPM 093		
G	DPM 278	DPM 278	DPM 278	DPM 278	DPM 278	DPM 278		
Н	DPM 128	DPM 129	DPM 130	DPM 131	DPM 132	DPM 133		

#### 6.7 AUTOCLAVING

DiamondAPEX pipettors are completely autoclavable at 1 bar pressure and 121°C for 20 minutes exposure time.

## **AUTOCLAVING INSTRUCTIONS**

- Keep the volume adjustment in the unlocked position.
- · Do not disassemble the pipettor for autoclaving.
- After autoclaving, allow the pipettor to completely cool and fully dry for a minimum of four (4) hours.

If the pipettor is autoclaved frequently, the piston and springs should be greased with the supplied lubricant to maintain smooth movement.

## 7 TROUBLE SHOOTING GUIDE

Problem area	Possible Cause	Solution
	Worn o-ring or seal	Replace worn parts
Pipettor is leaking	Foreign particles between tip and tip cone	Clean tip cone, attach new tip
	Foreign particles between piston and seal	Clean seal and piston
	Worn o-ring or seal	Replace worn parts
Pipettor does not	Tip cone is loose	Tighten tip cone
aspirate the solution	Piston is damaged (Chemically)	Return pipettor to authorised distributor
	Damaged tip cone	Replace the tip cone
	Improper assembly	See "Maintenance" section
Pipettor is	Tip cone is loose	Tighten tip cone
inaccurate	Tip incorrectly attached	Attach firmly
	Calibration altered	Recalibrate according to instructions
Inaccurate dispensing with certain liquids	Calibration not suitable for particular liquid	Recalibrate with the liquid in question

