

# GCC-E50 / Enhanced Clinical Centrifuge

## Quick Start Guide


### What's in the box:

- Quick Start Guide
- Centrifuge
- Power Cables (3 each)
- Rotor
- Wrench
- Door Key (small screwdriver)



### UNPACKING THE UNIT

Please be sure to remove any packing materials used during initial shipment:



1) Open the lid: Use the included Door Key to press the button on the right side of the unit.



2) Remove foam pads at the hinge latch and from around the rotor.

### LOADING INSTRUCTIONS

**IMPORTANT:** The rotor should always be loaded in a balanced manner to avoid personal injury and to ensure centrifuge longevity.

- When any tube is loaded, a tube of the same type, volume and weight must be added in the opposite rotor location.
- Never intentionally run the centrifuge in an unbalanced condition even if operation is not affected.

*Continued...*

# GCC-E50 / Enhanced Clinical Centrifuge


## Quick Start Guide

### OPERATION

**IMPORTANT:** Plug unit into appropriate main supply based on the rated voltage.

- 1 Select** ..... Press to choose the parameter you wish to modify.
- 2 Pulse** ..... Press to accelerate and hold at the set speed. The centrifuge stops immediately once this button is released.
- 3 Door** ..... Press to open the lid. This button is not enabled while the centrifuge is running.
- 4 Start / Stop** ..... Press to start the centrifuge. Press to stop if centrifuge is already running.
- 5 Adjustment Knob** Rotate clockwise to increase parameter values. Rotate counterclockwise to decrease parameter values. Press to switch between speed and RCF.

### CONTROL PANEL



Rotor Type	Tube Type	Tubes Per Rotor	Dimension	Maximum Speed (rpm)	Radius (cm)	Maximum RCF (xg)
CC-R2 (6x50mL & 6x15mL)	50mL Conical	6	29.5x100	4000	11	1800
	50mL Conical	6	17x120	4000	11	1800

For the full User Manual, please visit [www.globescientific.com](http://www.globescientific.com)



#### Customer Service, Sales and Technical Support:

Phone: 1 (800) 394-4562 • 1 (201) 599-1400

Fax: 1 (201) 599-1406 • E-mail: [mail@globescientific.com](mailto:mail@globescientific.com)