


Premium Microscope Slides

Globe Scientific's premium quality slides are quickly replacing overpriced, name-brand slides in markets throughout the world. From the smooth finish of our chemically etched frosted slides, to our unique safety corners, Globe slides will exceed your expectations.


Produced from premium glass, Globe slides are packaged in boxes of 72 and feature a moisture absorbent sheet in each box to prevent sticking. They are pre-cleaned and ready to use.



90° Corners
Traditional design



Safety Corners
Slightly trimmed to prevent finger and hand cuts



Clipped Corners
No sharp corners with 45° beveled edges. Perfect for blood smears.



Plain Slides

- Clear for full use of slide
- Dimensions: 25mm x 75mm x 1mm
- Thickness: 1 – 1.2mm
- 72/box, 20 boxes/case (10 Gross)


| Item | Description | Unit |
|------|--------------------------------|-------------|
| 1301 | Ground Edges, 90° Corners | 1440 slides |
| 1321 | Ground Edges, Safety Corners | 1440 slides |
| 1331 | Beveled Edges, Clipped Corners | 1440 slides |



Frosted Slides

- Chemically etched on one end, one side
- Smooth finish won't wear out markers
- Dimensions: 25mm x 75mm x 1mm
- Thickness: 1 – 1.2mm
- 72/box, 20 boxes/case (10 Gross)

| Item | Description | Unit |
|------|--------------------------------|-------------|
| 1304 | Ground Edges, 90° Corners | 1440 slides |
| 1324 | Ground Edges, Safety Corners | 1440 slides |
| 1334 | Beveled Edges, Clipped Corners | 1440 slides |



Two-Sided Frosted Slides

- Chemically etched on one end, two sides
- Smooth finish won't wear out markers
- Dimensions: 25mm x 75mm x 1mm
- Thickness: 1 – 1.2mm
- 72/box, 20 boxes/case (10 Gross)

| Item | Description | Unit |
|------|--------------------------------|-------------|
| 1308 | Ground Edges, 90° Corners | 1440 slides |
| 1328 | Ground Edges, Safety Corners | 1440 slides |
| 1338 | Beveled Edges, Clipped Corners | 1440 slides |

Frosted Color Slides


- Color is applied to one side of one end
- Dimensions: 25mm x 75mm x 1mm
- Thickness: 1 – 1.2mm
- 72/box, 20 boxes/case (10 Gross)

Ground Edge with 90° Corners

| Item | Color | Unit |
|-------|--------|-------------|
| 1304W | White | 1440 slides |
| 1304B | Blue | 1440 slides |
| 1304G | Green | 1440 slides |
| 1304N | Orange | 1440 slides |
| 1304P | Pink | 1440 slides |
| 1304T | Tan | 1440 slides |
| 1304Y | Yellow | 1440 slides |

Ground Edge with Safety Corners

| Item | Color | Unit |
|-------|--------|-------------|
| 1324W | White | 1440 Slides |
| 1324A | Aqua | 1440 Slides |
| 1324B | Blue | 1440 Slides |
| 1324D | Gold | 1440 Slides |
| 1324G | Green | 1440 Slides |
| 1324L | Lilac | 1440 Slides |
| 1324N | Orange | 1440 Slides |
| 1324P | Pink | 1440 Slides |
| 1324T | Tan | 1440 Slides |
| 1324Y | Yellow | 1440 Slides |



Beveled Edge with Clipped Corners


| Item | Color | Unit |
|-------|--------|-------------|
| 1334W | White | 1440 Slides |
| 1334B | Blue | 1440 Slides |
| 1334G | Green | 1440 Slides |
| 1334P | Pink | 1440 Slides |
| 1334Y | Yellow | 1440 Slides |



See reverse side for adhesive slides, slides with cavity wells and slide cover glass ▶

Globe Scientific also offers a complete line of slide mailers, storage boxes and staining accessories. Please call for details:
1(800) 394-4562

Premium Microscope Slides and Cover Glass



Adhesive Slides


- Tissue sections and cells adhere without adhesives
- Ideal for frozen sections and cytological preparations
- Integrated "+" symbol for easy identification
- Dimensions: 25mm x 75mm x 1mm
- Thickness: 1 – 1.2 mm
- 72/box, 20 boxes/case (10 Gross)

Ground Edge with 90° Corners

| Item | Color | Unit |
|-------|--------|-------------|
| 1358W | White | 1440 slides |
| 1358A | Aqua | 1440 slides |
| 1358B | Blue | 1440 slides |
| 1358G | Green | 1440 slides |
| 1358L | Lilac | 1440 slides |
| 1358N | Orange | 1440 slides |
| 1358P | Pink | 1440 slides |
| 1358T | Tan | 1440 slides |
| 1358Y | Yellow | 1440 slides |

Beveled Edge with Clipped Corners

| Item | Color | Unit |
|-------|--------|-------------|
| 1354W | White | 1440 Slides |
| 1354P | Pink | 1440 Slides |
| 1354Y | Yellow | 1440 Slides |



Slides with Cavity Wells

- Available in single or double cavity wells
- Well slides are ideal for holding liquids and cultures for microscope examination
- Made from selected 1.2mm thick microscope slides with 90° ground edges
- Spherically polished cavities are 13mm diameter with 0.5-0.6mm depth
- Dimensions: 25mm x 75mm x 1mm

| Item | Color | Unit |
|---------|--------------------------------------|------|
| 1341 | Single-cavity, 72/box, 20 boxes/case | 1440 |
| 1341-72 | Single-cavity, 72/box | 72 |
| 1344 | Double-cavity, 72/box, 20 boxes/case | 1440 |
| 1344-72 | Double-cavity, 72/box, 20 boxes/case | 72 |



Cover Glass Available in three thicknesses:

- No. 1: 0.14mm, +/- 0.015mm
- No. 1.5: 0.175mm, +/- 0.015mm
- No. 2: 0.21mm, +/- 0.02mm



Cover Glass

- Packaged 1 oz per box, 10 boxes per carton (10 oz)
- Each flip-top box contains a desiccant and is vacuum sealed in an aluminum wrap to prevent sticking

| Item | Description | Unit | Approximate Pcs per 1 oz Box |
|---------|--------------------|-------|------------------------------|
| 1401-10 | 18 x 18mm, No. 1 | 10 oz | 248 pieces |
| 1401-20 | 18 x 18mm, No. 2 | 10 oz | 168 pieces |
| 1402-10 | 20 x 20mm, No. 1 | 10 oz | 192 pieces |
| 1402-20 | 20 x 20mm, No. 2 | 10 oz | 154 pieces |
| 1404-10 | 22 x 22mm, No. 1 | 10 oz | 168 pieces |
| 1404-15 | 22 x 22mm, No. 1.5 | 10 oz | 148 pieces |
| 1404-20 | 22 x 22mm, No. 2 | 10 oz | 136 pieces |
| 1407-10 | 22 x 32mm, No. 1 | 10 oz | 114 pieces |
| 1407-20 | 22 x 32mm, No. 2 | 10 oz | 84 pieces |
| 1411-10 | 22 x 40mm, No. 1 | 10 oz | 86 pieces |
| 1411-20 | 22 x 40mm, No. 2 | 10 oz | 75 pieces |
| 1414-10 | 22 x 50mm, No. 1 | 10 oz | 70 pieces |
| 1414-20 | 22 x 50mm, No. 2 | 10 oz | 63 pieces |
| 1418-10 | 22 x 60mm, No. 1 | 10 oz | 60 pieces |
| 1418-20 | 22 x 60mm, No. 2 | 10 oz | 50 pieces |
| 1405-10 | 24 x 24mm, No. 1 | 10 oz | 141 pieces |
| 1405-20 | 24 x 24mm, No. 2 | 10 oz | 102 pieces |
| 1408-10 | 24 x 32mm, No. 1 | 10 oz | 104 pieces |
| 1408-20 | 24 x 32mm, No. 2 | 10 oz | 78 pieces |
| 1412-10 | 24 x 40mm, No. 1 | 10 oz | 82 pieces |
| 1412-20 | 24 x 40mm, No. 2 | 10 oz | 71 pieces |
| 1415-10 | 24 x 50mm, No. 1 | 10 oz | 65 pieces |
| 1415-15 | 24 x 50mm, No. 1.5 | 10 oz | 56 pieces |
| 1415-20 | 24 x 50mm, No. 2 | 10 oz | 56 pieces |
| 1419-10 | 24 x 60mm, No. 1 | 10 oz | 55 pieces |
| 1419-20 | 24 x 60mm, No. 2 | 10 oz | 46 pieces |



www.globescientific.com

Visit our website for the latest products
Customer service, sales and technical support

Phone 1 (800) 394-4562
1 (201) 599-1400

Fax 1 (201) 599-1406
E-mail mail@globescientific.com