

NEW Product Announcement

Carboys



Carboys – the tamper evident, sure seal caps on these carboys follow DIN recommendations

Carboys

Store, mix and dispense bulk quantities of liquids, solvents and cleaning solutions. Globe Scientific carboys are available in four sizes and feature a standardized industrial packaging design. They can be used in industrial packaging or for everyday laboratory practices. Produced from high density polyethylene (HDPE), they feature comfortable handles and tamper evident, secure seal caps that follow DIN recommendations.

The 5, 10 and 20 liter sizes conform to the UN regulation on transportation of dangerous liquids.

All are approved for use with food products.

Technical Information:

Item #	Volume	DIN	Dimensions
601780	2 Liter	GL 45	175 x 155 x 115mm
601781	5 Liter	GL 51	250 x 195 x 150mm
601782	10 Liter	GL 51	320 x 225 x 195mm
601783	20 Liter	GL 51	370 x 295 x 245mm

Ordering Information:

Item #	Description	Packaging Options
601780	2 Liter (.5 Gallon), HDPE	Each & 48/case
601781	5 Liter (1.3 Gallon), HDPE	Each & 30/case
601782	10 Liter (2.6 Gallon), HDPE	Each & 24/case
601783	20 Liter (5.3 Gallon), HDPE	Each & 6/case



Carboys with spigot on next page . . .

Carboys with Spigot



Carboys with Spigot



Carboys with Spigot

These heavy duty carboys are designed for lifetime service and can be used in sub-zero temperatures. They feature very thick walls but still retain translucency. The 3/4" BSP spigot (supplied) pairs with a 12mm deep thread in the carboy for total leak-proof performance. Note: Apply a few drops of Vaseline to the threading before assembly.

All are approved for use with food products.

Material: Heavy-duty polyethylene (HDPE)

Technical Information:

Item #	Neck I.D.	Neck O.D.	Bottle O.D.	Height
601660	52.5mm	62.5mm	170mm	345mm
601662	52.5mm	62.5mm	210mm	425mm
601664	79.5mm	95.5mm	280mm	565mm
601666	79.5mm	95.5mm	350mm	700mm

Ordering Information:

Item #	Description	Unit
601660	5 Liter (1.3 Gallon), HDPE	Each
601662	10 Liter (2.6 Gallon), HDPE	Each
601664	25 Liter (6.6 Gallon), HDPE	Each
601666	50 Liter (13.2 Gallon), HDPE	Each

globe
SCIENTIFIC INC.

www.globescientific.com

Globe Scientific Inc.

610 Winters Avenue • Paramus, NJ 07652

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com