

3mL Sample Cup, Polystyrene (PS)

Item# 110911



Description

These Globe Scientific multi-purpose sample cups are designed for use on the most popular instruments in the industry. The sample cups are made to exact specifications and provide excellent optical properties.

Features

- 3mL capacity
- Excellent clarity

Technical Specifications

Autoclavable: No
 Brittleness temperature: 20°C
 Maximum temperature: 90°C
 Sterility: Non-sterile

Products

Item#	Description	Unit
110911	3mL sample cup, PS	1000

