

Technical Data Sheet

Rev. 15348

Page 1 of 1

3mL Sample Cup, Polystyrene (PS)

Item# 110913



Products

| Item# | Description | Unit |
|--------|--------------------|------|
| 110913 | 3mL sample cup, PS | 1000 |

Description

These Globe Scientific sample cups are designed for use on the Tosoh® 360 analyzer. They are made to exact specifications and provide excellent optical clarity.


Features


- 3mL capacity
- Excellent clarity

Technical Specifications

Autoclavable: No
Brittleness temperature: 20°C
Maximum temperature: 90°C

Sterility: Non-sterile





Note: All instruments and names of manufacturers mentioned herein are the trademarks of their respective companies and are used for reference purposes only. Globe Scientific Inc. expressly disclaims and affiliation to them.



Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406 E-mail mail@globescientific.com Globe Scientific Inc.

400 Corporate Drive • Mahwah NJ 07430

www.globescientific.com

© Copyright 2015 Globe Scientific Inc. All rights reserved.