

12 x 57mm (3mL Capacity) Flat Bottom Polystyrene (PS) Tube

Item# 114010



Description

These Globe Scientific flat bottom, self-standing tubes are made from clear polystyrene (PS). The 3 mL tubes are graduated with rim.

Features

- 3mL capacity
- 12 x 57mm
- Polystyrene (PS)
- Graduated
- With rim
- Self-standing

Technical Specifications

Spin Rate: 1500 RPM
 Autoclavable: No
 Brittleness temperature: 20°C
 Maximum temperature: 90°C
 Sterility: Non-sterile

Products

Item#	Description	Unit
114010	3mL Flat bottom tube, PS	1000

