

Technical Data Sheet

Rev. 17006 Page 1 of 1

2mL Graduated Microcentrifuge Tubes with Attached Snap Cap, Polypropylene (PP)

Item# 111568



Description

Graduated microcentrifuge tube, 1000/pack. 2mL capacity. Attached leak-resistant snap cap. Molded graduations. These clear polypropylene tubes are resistant to most chemical reagents and can be autoclaved at 121°C (with cap open).

Features


- Polypropylene (PP)
- Attached snap cap, pierceable with writing area
- Molded graduations
- Frosted writing area on tube
- Non-sterile

Products

Item#	Description	Packaging
111568	2mL, PP, clear	1000 per pack

Technical Specifications

Material: Polypropylene (PP) Autoclavable: Yes. 121°C for 20 mintues (with cap open) Temperature range: --80°C to +121°C Spin rate: 16,400 rpm or 25,000 x g





 Customer service, sales and technical support:

 Phone
 1 (800) 394-4562 • 1 (201) 599-1400

 Fax
 1 (201) 599-1406

 E-mail
 mail@globescientific.com

Globe Scientific Inc. 400 Corporate Drive • Mahwah NJ 07430

www.globescientific.com © Copyright 2017 Globe Scientific Inc. All rights reserved.