

Exact Volume Transfer Pipet, 100 μ L (0.10mL)

Item# 139114

Technical Specifications

Material: Low Density Polyethylene (LDPE)
Length: 84mm (+/- 2mm)
Capacity: 100 μ L
Bulb Draw: 100 μ L
Sterility: Non-Sterile
Autoclavable: No
Temperature Range: -50°C to +75°C
Tolerance: +/- 15 μ L

Chemical Resistance

- Excellent resistance (no attack) to dilute and concentrated acids, alcohols, bases and esters
- Good resistance (minor attack) to aldehydes, ketones and vegetable oils
- Limited resistance (moderate attack suitable for short-term use only) to aliphatic and aromatic hydrocarbons, mineral oils, and oxidizing agents
- Poor resistance, and not recommended for use with halogenated hydrocarbons.

Quality Control

- This product is lot inspected and tested to ensure the highest level of quality.

Recommended Storage Conditions

- Store in cool dry place
- Avoid temperature variations and humidity

