

Wall Mount Drying/Draining Rack

Globe Scientific offers this sturdy 72 place laboratory drying rack. The rack is perfect for draining and drying all types of labware. It comes with a mounting kit to hang on the wall. The back plate of the support is made from a single mold of 4mm rigid high impact polystyrene (HIPS). The peg sockets have closed ends eliminating leakage and consequent biohazards. The seamless drip channel collects waste which drains via the supplied drain tube.

The 72 (push in) pegs can accommodate items with a neck bore in excess of 15mm. They are easily removable for cleaning or to accommodate odd shaped items. Also included are 11 small pegs to accommodate small diameter items, such as test tubes. This durable rack will not rust or stain.

Features:

- Heavy duty construction
- 72 interchangeable pegs
- Drip channel with draining hose
- Mounting kit and template for hanging
- Hooks for joining together additional racks

Dimensions: 24.75" H x 17.75" W x 4.5" D



Item#	Size	Unit
600213	Drying rack, 72-place	Each
601213	Additional small pegs, 11/pack	11