

# Borosilicate Glass Culture Tubes

Globe Scientific offers this outstanding line of high quality **borosilicate glass culture tubes**. These tubes offer exceptional quality at great prices. They are neatly packaged in heavy duty boxes and shrink wrapped to keep them particle free. The boxes are compact and designed to fit into most test tube dispensers.

### Features:

- Tubes can be used for tissue culture and blood banks.
- Tubes are chemically inert.
- Made from superior quality borosilicate glass to reduce pH changes and contaminants that are potentially leached from tubes made from soda lime glass.
- Tubes feature well-rounded bottoms and smooth, fire-polished rims.
- Dimensional specifications are rigidly controlled for reliable, uniform results.
- Produced in an ISO 9001:2008 certified factory.

Item#	Size (mm)	Overflow Capacity	Packaging	Unit
1503	10 x 75	4mL	250/box, 4 boxes/case	1000
1505	12 x 75	6mL	250/box, 4 boxes/case	1000
1510	13 x 100	10mL	250/box, 4 boxes/case	1000
1512	16 x 100	14mL	250/box, 4 boxes/case	1000
1515	16 x 125	19mL	250/box, 4 boxes/case	1000
1517	16 x 150	23mL	250/box, 4 boxes/case	1000
1519	18 x 150	29mL	125/box, 4 boxes/case	500
1522	20 x 150	36mL	125/box, 4 boxes/case	500



Globe borosilicate glass culture tubes are available in 8 sizes.



Boxes are compact and designed to fit into most test tube dispensers.

### Dispensers for Glass Culture Tubes

Made from heavy duty white painted metal. Front plexi-glass cover pulls up to dispense the tubes.

Item#	Description
1590	for 10 x 75mm and 12 x 75mm tubes
1592	for 13 x 100mm tubes



[www.globescientific.com](http://www.globescientific.com)

Globe Scientific Inc.

610 Winters Avenue

Paramus, NJ 07652

**Customer service, sales and technical support:**

Phone 1 (800) 394-4562

1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail [mail@globescientific.com](mailto:mail@globescientific.com)