

Mini Tube Rocker



MINI TUBE ROCKER

Globe Scientific's mini tube rocker is designed to provide thorough agitation of tube sizes ranging from 1.5mL to 50mL. A grooved rubber mat securely holds the tubes in place.

Small in size (only 5 x 13 inches), the instrument can fit on any lab bench and can also be placed inside a cold room, refrigerator or incubator.


Item #	Description
527	Mini Tube Rocker, 12" x 4.5" platform, 115V
527-230	Mini Tube Rocker, 12" x 4.5" platform, 230V

SPECIFICATIONS:

Platform Size:	12" x 4.5" (30cm x 11.4cm)
Speed - 115V:	Fixed at 24 rpm
Speed - 230V:	Fixed at 20 rpm
Tilt Angle:	30 Degrees
Work Surface:	Grooved rubber mat
Operating Temp. Range:	+4°C to 65°C
Dimensions:	13.3" x 5.3" x 4"
Weight:	4 lbs (2 kgs)
Electrical:	115V or 230V, 50-60Hz, 50W
Warranty:	2 Years

- Small footprint - fits on any lab bench
- Ideal for all tube sizes 1.5mL to 50mL
- Optimized speed and tilt angle for thorough mixing in sample tubes
- Cold room and incubator safe
- No adapters or mounting hardware required

Ideal for a variety of tube sizes:



globe[®]
SCIENTIFIC INC.

www.globescientific.com

Globe Scientific Inc.

610 Winters Avenue • Paramus, New Jersey 07652 • USA

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com