

Sample Tubes for Storage and Transportation

Globe Scientific's sample tubes have been designed for the storage and transport of biological material. They are produced from non-toxic polypropylene (PP) and offer a high degree of strength and clarity. The tubes feature external threads which

provide a smooth, uniform inner surface.

- Tubes withstand a temperature range of -196 °C to +121 °C
- Round bottom tubes can be centrifuged up to 17,000 x g
- Autoclavable (121 °C) in an upright position with caps loosened

Non-Graduated

Caps sold separately. See reverse for details.

(shown actual size)



Order Number	6030	6031	6032	6033	6034	6035	6038	6039
Volume	1.2 mL	2.0 mL	2.0 mL	3.0 mL	4.0 mL	4.0 mL	5.0 mL	5.0 mL
Height	40.9 mm	44.7 mm	47.4 mm	69.4 mm	72.5 mm	74 mm	88 mm	89 mm
Style	Self standing	Round bottom	Self standing	Self standing	Round bottom	Self standing	Round bottom	Self standing

Graduated with White Writing Area

Caps sold separately. See reverse for details.

(shown actual size)



Order Number	6050	6051	6052	6053	6054	6055	6058	6059
Volume	1.2 mL	2.0 mL	2.0 mL	3.0 mL	4.0 mL	4.0 mL	5.0 mL	5.0 mL
Height	40.9 mm	44.7 mm	47.4 mm	69.4 mm	72.5 mm	74 mm	88 mm	89 mm
Style	Self standing	Round bottom	Self standing	Self standing	Round bottom	Self standing	Round bottom	Self standing



Screw Caps


The screw caps for Globe Scientific's sample tubes are available either with or without an o-ring.

Screw Cap without O-Ring

Order Number	Color
6136	Natural ○
6136B	Blue ●
6136A	Gray ●
6136G	Green ●
6136L	Lilac ●
6136N	Orange ●
6136R	Red ●
6136W	White ○
6136Y	Yellow ●

Screw Cap with O-Ring

Order Number	Color
6137	Natural ○
6137B	Blue ●
6137A	Gray ●
6137G	Green ●
6137L	Lilac ●
6137N	Orange ●
6137R	Red ●
6137W	White ○
6137Y	Yellow ●



Workstation Racks

Workstation securely holds 40 vials each. The Globe Scientific self-standing vials interlock into the rack with a simple twist of the hand.

- 40-place format (10 x 4)
- Made from polycarbonate with anti-slip rubber feet
- Workstation racks are designed to securely lock vials into place which allows vials to be opened or securely closed with one hand
- Withstand temperature range of -196°C to +121°C
- Autoclavable at 121°C, 15 psig for 20 minutes

Dimensions: 21.5 x 10.0 x 2.50 cm

Packaging: 5 units/box

Order Number	Color
3049R	Red
3049B	Blue
3049G	Green



www.globescientific.com

Visit our website for the latest products
Customer service, sales and technical support

Phone 1 (800) 394-4562

1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com

Globe Scientific Inc.
 610 Winters Ave.
 Paramus, NJ 07652