

Graduated Transfer Pipet, 7mL

Item# 135030

Technical Specifications

Material:	Low Density Polyethylene (LDPE)
Length:	160mm (+/- 0.5mm)
Capacity:	7mL
Drops/mL:	20
Bulb Draw:	3.2mL
Graduation Increments:	0.5mL
Graduated to:	3mL
Sterility:	Non-Sterile
Autoclavable:	No
Temperature Range:	-50°C to +75°C

Chemical Resistance

- Excellent resistance (no attack) to dilute and concentrated acids, alcohols, bases and esters
- Good resistance (minor attack) to aldehydes, ketones and vegetable oils
- Limited resistance (moderate attack suitable for short-term use only) to aliphatic and aromatic hydrocarbons, mineral oils, and oxidizing agents
- Poor resistance, and not recommended for use with Halogenated hydrocarbons.

Quality Control

- This product is lot inspected and tested to ensure the highest level of quality.

Recommended Storage Conditions

- Store in cool dry place
- Avoid temperature variations and humidity

