

The Best Microscope Slide Storage Cabinet Ever!



Image shows base and 3 cabinet units.


Dimensions:

A-Cabinet Unit: 402mm wide x 477mm deep x 126mm high
(15.75" wide x 18.75" deep x 5" high)

B-Base: 402mm wide x 477mm deep x 129mm high
(15.75" wide x 18.75" deep x 5" high)



Brass pull knobs and card holders on each dual slotted drawer. Included foam inserts keep slides neatly standing.



Interlocking system makes it possible to safely store up to 45,000 slides on a single stack.

This sturdy, heavy-duty metal slide storage cabinet is the ideal solution for the safe, long-term storage of slide specimens. Each cabinet unit holds up to 4,500 slides. Up to 10 units can be secured together on a single base for a total of 45,000 slides per stack!

- Each cabinet unit contains 6 dual slotted drawers
- Each drawer holds up to 750 slides (4,500 per 6 drawer cabinet unit)
- Expandable: safely stack up to 10 cabinet units on a single base (total of up to 45,000 slides per stack)
- One base is required for each new stack of up to 10 cabinets
- Interlocking system locks the stacked cabinet units together
- Each drawer has brass pull knob and identification card holder — identification cards included
- Included foam blocks keep slides upright in drawers
- Fits all standard size microscope slides: 25 x 75mm and 1" x 3"
- Chemical resistant plastic powder coating
- Smooth sliding drawers
- Available in green, tan and white

Item #	Description	Unit
513500G	Slide Storage Cabinet, Green	Each
513503G	Slide Storage Base, Green	Each
513500T	Slide Storage Cabinet, Tan	Each
513503T	Slide Storage Base, Tan	Each
513500W	Slide Storage Cabinet, White	Each
513503W	Slide Storage Base, White	Each

globe
SCIENTIFIC INC.

Customer service, sales and technical support:

Phone: 1 (800) 394-4562 • 1 (201) 599-1400 • Fax: 1 (201) 599-1406
E-mail: mail@globescientific.com • www.globescientific.com