

Diamond[®] Jr. Pocket-Size Pipettor

Diamond[®] Jr. pipettors are pocket-sized versions of our high-precision Diamond[®] pipettors. Designed for use in diagnostic kits, water testing, science classrooms and all laboratories, these pipettors provide an economical and accurate method for pipetting fixed volumes.

- Fits easily in a lab coat pocket
- Two-stage plunger for last-drop dispensing
- Ideal for diagnostic kits where an economical pipettor is required
- Completely autoclavable
- Color-coded for easy identification
- Manufactured in an ISO 9001:2015 facility

Item #	Volume (µl)	Color	Unit
3322-5	5	Red	Each
3322-10	10	Pink	Each
3322-20	20	Green	Each
3322-25	25	Green	Each
3322-35	35	Blue	Each
3322-50	50	Purple	Each
3322-100	100	Orange	Each
3322-200	200 <i>New Volume</i>	Yellow	Each



3322-200

NEW VOLUME



www.globescientific.com

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com