

**FULL TURN Screw Cap**

# 4 oz Urine Collection Containers

**ALL NEW!**



**Designed for on-site collection and testing.**

These 4 ounce urine collection containers are ideal for use in physician's offices and hospital settings. The containers feature a full turn screw cap for better leak resistance. They have molded graduation marks and a frosted writing area. The individually wrapped sterile containers also come complete with an affixed patient ID label.

- Graduated to 120mL
- Frosted writing area
- Clarified polypropylene (PP)

**Sterile & Aseptic Containers:**

- Individually wrapped
- Trilingual patient ID label (English, Spanish & French)

Material: Container = Polypropylene (PP) • Cap = Polyethylene (PE)  
Dimensions: Container: 57.6 x 72.0mm • Cap: 62mm O.D.

Item #	Description	Unit
5901G	Full turn green attached screw cap, ID label, individually wrapped, STERILE	250
5901R	Full turn red attached screw cap, ID label, individually wrapped, STERILE	250
5902G	Full turn green attached screw cap, ID label, individually wrapped, CLEANROOM ISO8 ASEPTIC	250
5902R	Full turn red attached screw cap, ID label, individually wrapped, CLEANROOM ISO8 ASEPTIC	250
5905G	Full turn green separate screw cap, non-sterile	500
5905R	Full turn red separate screw cap, non-sterile	500

**diamond**ESSENTIALS  
*Quality. Value. No Compromises.*