

# GTR-2-FS / Fixed Speed Two-Tier Tube Rocker



*GTR-2-FS  
Fixed Speed Two-Tier Tube Rocker*

The GTR-2-FS fixed-speed tube rocker takes the proven design of the GTR-FS single tier model and adds another tier for even more capacity. With its see-saw motion, a fixed speed of 24rpm and an angle of +/-25° it provides optimal mixing. The ridged platform pad design allows the operator to place any tube from 1.5mL to 50mL on either tier of the unit, turn the power switch on and walk away.

**Features:**

- Adds an additional platform tier for double the capacity
- Simple plug & play operation
- Ridge pad design for easy, stable tube placement
- Small size uses minimal bench space
- Accommodates tubes from 1.5mL to 50mL

***When a single platform isn't enough,  
choose the space-saving design of the  
Globe Scientific GTR-2-FS!***

# Fixed Speed Two-Tier Tube Rocker



Accommodates (12) 50mL Tubes

Also ideal for a variety of additional tube sizes:

(24) 15mL Tubes

(48) 5mL Tubes

(24) Blood Collection Tubes (5mL, 7mL, 10mL or 15mL)

(27) Microcentrifuge Tubes (1.5mL or 2.0mL)

## GTR-2-FS Specifications

Working area (each tier)	252mm W x 170mm D (9.9in W x 6.7in D)
Maximum load	7lbs (3.2kg)
Maximum tilt angle	25°
Voltage	120v 60Hz or 240v 50Hz
Power	6W
Speed	24rpm fixed
Dimensions (WxDxH)	346x196x170mm / 13.6 x 7.7 x 7in
Weight	2.7kg / 6lbs
Operating temperature	5-50 °C
Permissible relative humidity	80%RH (non-condensing)
Warranty	2 Years

## ORDERING INFORMATION

Item #	Description	Unit
GTR-2-FS	Tube Rocker, double tier, 120v, 60Hz, US plug, 9.9" x 6.7" platforms, 25rpm	Each
GTR-2-FS-UK-EU	Tube Rocker, double tier, 240v, 50Hz, UK & EU plugs, 9.9" x 6.7" platforms, 25rpm	Each

If two platforms is more than you need, ask about the Globe Scientific GTR-FS Single Tier Rocker!



GTR-FS



[www.globescientific.com](http://www.globescientific.com)

### Customer Service, Sales and Technical Support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail [mail@globescientific.com](mailto:mail@globescientific.com)