

GTR-DVS / Variable Speed

Digital Tube Roller



a wide range of tubes, containers and bottles. The speed can be set digitally anywhere between 10rpm and 70rpm. In addition to a rolling motion, the rollers are mounted slightly off-center to rise and fall in a "see-saw" fashion. These two motions ensure your samples are thoroughly mixed for optimal homogeneity. The six rollers can accommodate up to 8.8lbs or ten 15ml/ 50ml centrifuge tubes.

Features:

- LCD display of speed and time
- Speed range of 10-70rpm
- Six roller design for gentle sample agitation and ease of cleaning
- Continuous or timed operation
- Rolling motion with see-saw movement
- View both speed and time concurrently
- Cold room and incubator compatible
- DC brushless motor is quiet and maintenance free

GTR-DVS Digital Tube Roller

For simple, accurate operation and thorough mixing of your tubes, bottles and containers, choose the Globe Scientific GTR-DVS!

The GTR-DVS can accurately set time and speed using the simple parameter adjustment knob and digital display



The six roller configuration can accommodate up to ten 15mL or ten 50mL tubes and is compatible with a variety of other tube sizes

Variable Speed Digital Tube Roller

Specifications	GTR-DVS	
Amplitude	24mm	
Motion	Rocking and rolling	
Max. Load Capacity	4kg (8.8lbs)	
Number of Rollers	6	
Roller Size (Length)	280mm (11in)	
Speed Range	10 -70rpm	
Speed Display	LCD	
Timer	Yes	
Timer Display	LCD	
Time Setting Range	1-1199min	
Operation Mode	Continuous or timed	
Voltage	100-240V, 50/60Hz	
Power	30W	
Dimensions (WxDxH)	260 x 450 x 120mm (10.2 x 17.7 x 4.7in)	
Weight	5.1kg (11.2lbs)	
Operating Temperature	5° - 40° C	
Permissible Relative Humidity	80% RH (non-condensing)	
Protection Class	IP21	
Warranty	2 Years	

ORDERING INFORMATION

Item#	Description	Unit
GTR-DVS	Variable Speed Digital Tube Roller, 6 Rollers, 120-240v, 50/60Hz	Each

www.globescientific.com

Customer Service, Sales and Technical Support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

1 (201) 599-1406 E-mail mail@globescientific.com

