

GTR-FS / Fixed Speed

Single Tier Tube Rocker



The GTR-FS fixed-speed tube rocker is the best choice for basic laboratory tube mixing applications. With its see-saw motion, a fixed speed of 24rpm and an angle of +/-25° it provides optimal mixing. The ridged platform pad design allows the operator to place any tube from 1.5mL to 50mL on the unit, turn the power switch on and walk away.

Features:

- Simple plug & play operation
- Ridge pad design for easy, stable tube placement
- Small size uses minimal bench space
- Accommodates tubes from 1.5mL to 50mL

GTR-FS Fixed Speed Single Tier Tube Rocker accommodates (6) 50mL tubes. See reverse side for additional tube size accommodations.

Fixed Speed Single Tier Tube Rocker


Also ideal for a variety of additional tube sizes:



Accommodates (6)
50mL Tubes



(12) 15mL Tubes



(24) 5mL Tubes



(12) Blood Collection
Tubes (5mL, 7mL, 10mL or 15mL)



(36) Microcentrifuge
Tubes (1.5mL or 2.0mL)

GTR-FS Specifications

Working area	252mm W x 170mm D (9.9in W x 6.7in D)
Maximum load	7lbs (3.2kg)
Maximum tilt angle	25°
Voltage	120v 60Hz or 240v 50Hz
Power	6W
Speed	24rpm fixed
Dimensions (WxDxH)	346x196x95mm / 13.6 x 7.7 x 3.7in
Weight	2.5kg / 5.5lbs
Operating temperature	5-50 °C
Permissible relative humidity	80%RH (non-condensing)
Warranty	2 Years

ORDERING INFORMATION

Item #	Description	Unit
GTR-FS	Tube Rocker, single tier, 120v, 60Hz, US plug, 9.9" x 6.7" platform, 25rpm	Each
GTR-FS-UK-EU	Tube Rocker, single tier, 240v, 50Hz, UK & EU plugs, 9.9" x 6.7" platform, 25rpm	Each

If a single platform isn't enough, ask about the space-saving design of Globe Scientific's GTR-2-FS Two-Tier Tube Rocker!



GTR-2-FS



www.globescientific.com

Customer Service, Sales and Technical Support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com