

GTRM-FS

VacuumTube Rockin' Roller

Quick Start Guide

What's in the box:

- Rockin' Roller
- Quick Start Guide

OPERATION

- Simply switch the main power to the rocker ON or OFF. When the unit is turned ON, it will roll at a fixed speed of 36rpm with a swing angle of $\pm 11^\circ$.



Power Switch

Located on right side of unit

IMPORTANT:

- The tube rocker/roller must not be overloaded or damage may occur. Please do not load the unit in excess of 6.6lbs (3kg).

www.globescientific.com

globe[®]
SCIENTIFIC INC.

Customer Service, Sales and Technical Support:

Phone: 1 (800) 394-4562 • 1 (201) 599-1400

Fax: 1 (201) 599-1406 • E-mail: mail@globescientific.com