

GVM-MINI / Fixed Speed Mini Vortex Mixer



*For basic vortexing in a small package,
choose the GVM-MINI!*

GVM-MINI Specifications

Voltage	100-240V,50/60Hz
Power	12W
Mixing motion	Orbital
Orbital diameter	4.8mm
Motor rating output	10W
Speed	3000rpm (fixed)
Run type	Touch operation
Dimensions[WxDxH]	133×133×80mm / 5.2 x 5.2 x 3.1in
Weight	0.6kg / 1.3lbs
Operating Temperature	5 to 40°C
Permissible relative humidity	80%RH (non-condensing)
Warranty	2 Years

When looking for a basic, fixed speed vortex mixer with a small footprint for your limited bench space, look no further than the Globe Scientific GVM-MINI! Despite its small size, the GVM-MINI has a wide, stable base to remain in place when in use. When you press your tube down on the top cup, the 3000rpm and small 4.8mm orbit create a 'vibratory' motion to quickly and effectively mix most tube sizes.

Features:

- Touch operation
- Fixed speed control of 3000rpm
- Used for mixing various tube sizes and formats
- Small footprint yet extremely stable
- Rubber suction feet ensure placement and mitigate vibration
- Spill resistant chassis

ORDERING INFORMATION

Item #	Description	Unit
GVM-MINI	Fixed Speed Mini Vortex Mixer, 120v, 60Hz, US Plug	Each
GVM-MINI-UK	Fixed Speed Mini Vortex Mixer, 240v, 50Hz, UK Plug	Each
GVM-MINI-EU	Fixed Speed Mini Vortex Mixer, 240v, 50Hz, EU Plug	Each