

Inoculation Loops



Loops feature integrated needle at reverse end



Packaged in resealable, tamper-evident bags

INOCULATION LOOPS - FLEXIBLE FEEL IN TAMPER-EVIDENT RESEALABLE BAGS

Globe Scientific's disposable inoculation loops and needles offer the most convenient and safest way of streaking culture media. These polystyrene (PS) inoculation loops are smooth and flexible to facilitate uniform streaking of culture media. The loops are packaged in tamper-evident resealable bags. All loops are sterilized by gamma irradiation. The 1µL and 10µL inoculation loops are available with the loop on one end and an integrated needle on the other end eliminating the need for a separate needle.

- Manufactured under ISO 9001:2008 and ISO 14001:2004 standards
- Gamma sterilized, sterility assurance level (SAL) 10^{-3}
- Certificate of calibration
- Tamper-evident resealable bags
- Integrated loop/needle design
- High impact polystyrene (HIPS)

Item #	Description	Unit
2871-01	1µL with needle, blue, PS, STERILE, individually wrapped	250
2871-10	1µL with needle, blue, PS, STERILE, 10/zipper seal bag (100 bags per case)	1000
2871-25	1µL with needle, blue, PS, STERILE, 25/zipper seal bag (40 bags per case)	1000
2875-01	10µL with needle, yellow, PS, STERILE, individually bagged	250
2875-10	10µL with needle, yellow, PS, STERILE, 10/zipper seal bag (100 bags per case)	1000
2875-25	10µL with needle, yellow, PS, STERILE, 25/zipper seal bag (40 bags per case)	1000



www.globescientific.com

Globe Scientific Inc.

400 Corporate Drive
Mahwah, New Jersey 07430 • USA

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com