

Lab Support Jacks

Rugged, heavy-duty aluminum alloy jacks provide stable height adjustment for general lab devices such as hot plates, water baths, glassware and more.

Features:

- High load capacity
- Stable, rigid construction
- Anodized aluminum with attractive satin finish
- Rust-proof, corrosion and chemical resistant for years of trouble-free performance
- Easy grip knurled adjustment knob
- Fully autoclavable at 121°C
- Easy to clean



Item #	Size	Min. to Max. Height	Load Capacity	Color	Unit
450000SM	80mm x 76mm (3.2" x 3.0")	48mm to 148mm (1.9" x 5.8")	30kg (66 lbs.)	Blue	Each
450000MD	140mm x 120mm (5.5" x 4.7")	58mm to 242mm (2.3" x 9.5")	50kg (110 lbs.)	Red	Each
450000LG	150mm x 121mm (5.9" x 4.8")	78mm to 278mm (3.1" x 10.9")	80kg (176 lbs.)	Yellow	Each