

NEW PRODUCT

NEW!



OIL CENTRIFUGE TUBES

Globe Scientific Oil Centrifuge Tubes facilitate testing of petroleum products for water and sediment. Available in both 6-inch "Short Cone" and 8-inch versions, these tubes match most equipment and test protocols. Premium 3.3 Borosilicate glass material insures strength, durability, and chemical resistance.


- 3.3 Borosilicate Glass Construction
- Graduated to 100mL in specified increments
- Blue Enamel Graduations
- White printed marking spot
- Re-usable
- Autoclavable

GRADUATION INTERVALS

Volume Range (mL)	Interval (mL)
0 to 0.5	0.05
0.5 to 2.0	0.1
2 to 3	0.2
3 to 5	0.5
5 to 10	1
10 to 25	5
25 to 50	25
50 to 100	50

ORDERING INFORMATION

Item #	Description	Length (mm)	Diameter (mm)	Unit(s)
8600106	Glass, Conical Bottom, 6 Inch, 100mL	165	45	Each & 10
8600108	Glass, Conical Bottom, 8 Inch, 100mL	200	36	Each & 10



Tubes shown actual size



www.globescientific.com

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com