



# Diamond RV-Pette-PRO™ Dispenser Syringe Tips

Precision manufacturing under strict quality controls ensures that RV-Pette PRO™ Dispenser Syringe Tips deliver excellent results when used with Diamond RV-Pette PRO™ Repeat Volume Pipettors and most popular brands of repeat volume pipettors. RV-Pette PRO™ tips are available in 12 sizes for dispensing volumes from 1µL to 5500µL and are available non-sterile and sterile.

## FEATURES:

**Variety** – Wide range of sizes match all sizes offered by other manufacturers and eliminates the need to change protocols or practice.

**Simplicity** – New colored adapter ring for automatic tip recognition on select Eppendorf® pipettors. Industry standard color coding of ring ensures easy volume identification.

**Safety** – Tips use the direct displacement method and are the safe, secure choice for handling a wide range of liquids including those that are viscous, infectious, and volatile.

**Reliability** – Each tip is tested during the automated manufacturing process.

**Confidence** – Sterile tips are individually wrapped and certified to be RNase, DNase, ATP, and pyrogen free for even the most critical applications

Non-Sterile	Sterile		
Item #	Item #	Size	Unit
3920	3920S	0.05mL	100
3921	3921S	0.1mL	100
3922	3922S	0.2mL	100
3923	3923S	0.5mL	100
3924	3924S	1.0mL	100
3925	3925S	1.25mL	100
3926	3926S	2.5mL	100
3927	3927S	5mL	100
3928	3928S	10mL	100
3929	3929S	12.5mL	100
3930	3930S	25mL*	25
3931	3931S	50mL*	25

\*Includes 1 tip adapter per box of 25.

Item #	Description	Unit
3930-A	Adapter for 25mL Tip	Each
3931-A	Adapter for 50mL Tip	Each

Item #	Description	Unit
3915	Dispenser Tip Starter Pack: (5) each tip size, non-sterile, (1) 25mL tip adapter, (1) 50mL tip adapter (1) dosage chart	Each



3920



3921



3922



3923



3924



3925




3926




3927




3928



3929



3930




3931

RV-Pette PRO™ Syringe Tip compatibility chart on reverse >>

# Diamond RV-Pette-PRO™ Dispenser Syringe Tips

## Compatibility Chart

The wide range of sizes match all sizes offered by other manufacturers and eliminates the need to change protocols or practice!



12  
tip sizes!

RV-Pette PRO dispenser tips feature automatic tip recognition on select Eppendorf® pipettors!

RV-Pette PRO  
Dispenser  
Syringe  
Tip Size

### PIPETTOR BRAND/MODEL:

	Diamond RV-Pette™	Diamond RV-Pette PRO™	Bel-Art Roxy-M™ Repeating	Biohit® Mechanical Stepper	BrandTech™ HandyStepper™ S	Capp R10™	Capp Rhythm™ Repeater	eLine® Light Dispenser	eLine® Pro Dispenser	Eppendorf® Repeater 4780	Eppendorf® Multipipette® plus	Fisherbrand™ Multipipette stream® /Xstream®	Fisherbrand™ HandyStep™ S	Gilson® Distriman®	Nichiryo NichiMate	Rainin® AutoRep S	Thermo® Fimppipette Stepper	VWR® Repeating Dispenser
0.05mL	YES	YES	✓	✓	–	✓	✓	–	–	✓	–	–	–	✓	–	–	–	✓
0.1mL	YES	YES	–	–	✓	✓	✓	✓	–	✓	✓	✓	–	–	✓	–	–	✓
0.2mL	YES	YES	–	–	–	✓	✓	✓	–	✓	✓	–	–	–	–	–	–	✓
0.5mL	YES	YES	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	–	✓	✓	✓	✓	✓
1.0mL	YES	YES	–	–	✓	✓	✓	✓	–	✓	✓	✓	–	–	✓	–	–	✓
1.25mL	YES	YES	✓	✓	✓	✓	✓	–	–	✓	–	–	✓	✓	✓	✓	✓	✓
2.5mL	YES	YES	–	✓	✓	✓	✓	✓	✓	✓	✓	✓	–	✓	✓	✓	✓	✓
5mL	YES	YES	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	–	✓	✓	✓	✓	✓
10mL	YES	YES	–	–	✓	✓	✓	✓	–	✓	✓	✓	–	–	✓	–	–	✓
12.5mL	YES	YES	✓	✓	✓	✓	✓	–	–	✓	–	–	✓	✓	✓	✓	✓	✓
25mL	YES	YES	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	–	✓	✓	✓	✓	✓
50mL	YES	YES	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	–	✓	✓	✓	✓	✓

All instruments and names of manufacturers mentioned herein are the trademarks of their respective companies and are used for reference purposes only. Globe Scientific Inc. expressly disclaims any affiliation to them.

Also Available:

Diamond RV-Pette-PRO  
Repeat Volume Pipettor

Call today or visit  
[www.globescientific.com](http://www.globescientific.com)  
for full product information.



**globe**  
SCIENTIFIC INC.

Customer service, sales and technical support:

Phone: 1 (800) 394-4562 • 1 (201) 599-1400 • Fax: 1 (201) 599-1406

E-mail: [mail@globescientific.com](mailto:mail@globescientific.com) • [www.globescientific.com](http://www.globescientific.com)

© Copyright 2019 Globe Scientific Inc. All rights reserved. rv-pettePRO-tips 19043