
Sample Tubes for Storage and Transportation



Globe Scientific's sample tubes have been designed for the storage and transport of biological material. They are produced from non-toxic polypropylene (PP) and offer a high degree of strength and clarity. The tubes feature external threads which

provide a smooth, uniform inner surface.

- Tubes withstand a temperature range of -196 °C to +121 °C
- Round bottom tubes can be centrifuged up to 17,000 x g
- Autoclavable (121 °C) in an upright position with caps loosened



Graduated with White Writing Area Caps sold separately. See reverse for details. (shown actual size)

| Order Number | 6050 | 6051 | 6052 | 6053 | 6054 | 6055 | 6058 | 6059 |
|--------------|---------------|--------------|---------------|---------------|--------------|---------------|--------------|---------------|
| Volume | 1.2 mL | 2.0 mL | 2.0 mL | 3.0 mL | 4.0 mL | 4.0 mL | 5.0 mL | 5.0 mL |
| Height | 40.9 mm | 44.7 mm | 47.4 mm | 69.4 mm | 72.5 mm | 74 mm | 88 mm | 89 mm |
| Style | Self standing | Round bottom | Self standing | Self standing | Round bottom | Self standing | Round bottom | Self standing |



Screw Caps

The screw caps for Globe Scientific's sample tubes are available either with or without an o-ring.

Screw Cap without 0-Ring

| Packaging | 1000/bag |
|--------------|------------------|
| Order Number | Color |
| 6136 | Natural O |
| 6136B | Blue |
| 6136A | Gray |
| 6136G | Green |
| 6136L | Lilac |
| 6136N | Orange 🛑 |
| 6136R | Red • |
| 6136W | White \bigcirc |
| 6136Y | Yellow — |

Screw Cap with 0-Ring

| Packaging | 1000/bag |
|--------------|------------------|
| Order Number | Color |
| 6137 | Natural O |
| 6137B | Blue |
| 6137A | Gray |
| 6137G | Green |
| 6137L | Lilac |
| 6137N | Orange 🛑 |
| 6137R | Red 🛑 |
| 6137W | White \bigcirc |
| 6137Y | Yellow 🛑 |



Workstation Racks

Workstation securely holds 40 vials each. The Globe Scientific self-standing vials interlock into the rack with a simple twist of the hand.

- 40-place format (10 x 4)
- Made from polycarbonate with anti-slip rubber feet
- Workstation racks are designed to securely lock vials into place which allows vials to be opened or securely closed with one hand
- Withstand temperature range of -196°C to +121°C
- Autoclavable at 121°C, 15 psig for 20 minutes

Dimensions: 21.5 x 10.0 x 2.50 cm

Packaging: 5 units/box

| Order Number | Color | | |
|--------------|-------|--|--|
| 3049R | Red | | |
| 3049B | Blue | | |
| 3049G | Green | | |



www.globescientific.com

Visit our website for the latest products Customer service, sales and technical support

Phone 1 (800) 394-4562 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com

Globe Scientific Inc. 610 Winters Ave. Paramus, NJ 07652