

Diamond[®] SeroFlow[™] and Diamond[®] SeroFlow[™] Plus Serological Pipette Controllers

Choose from two models to fit your lab's needs and budget.

diamond SeroFlow[™]

The **Diamond SeroFlow** has a comfortable, ergonomic design for hours of continuous use. An easy to adjust analog dial sets the coarse aspiration speed to work best with your method and use. Precise control is then achieved by pressing lightly or heavily on the triggers. An easy-to-replace silicone pipette holder and 0.45µm filter make maintenance a breeze.



3385 –
Diamond[®] SeroFlow[™]

diamond SeroFlow[™] plus

The **Diamond SeroFlow Plus** takes the solid foundation of the standard model and adds an even more powerful motor for faster aspiration speeds as well as a 6 speed digital coarse setting for the ultimate in consistency for varying solution viscosities, as well as an LED display showing the coarse speed setting and battery charge status.



3387 –
Diamond[®] SeroFlow[™] Plus



LED display on the Diamond SeroFlow Plus shows the coarse speed setting and battery charge status.



Diamond® SeroFlow™ and Diamond® SeroFlow™ Plus Serological Pipette Controllers



Features:

- Simple plug & play operation
- Ergonomic, lightweight designs to help reduce fatigue
- Powerful motors for fast liquid draw
- Precise control of coarse and fine speeds allows for both speed and accuracy
- Accommodates all size serological pipettes, glass and plastic
- 0.45µm hydrophobic filter prevents contamination
- Long life Lithium Ion battery for hours of use and no 'memory'
- Can use while charging, no waiting
- Autoclavable nose cone ensures contamination-free pipetting

Specifications	SeroFlow™	SeroFlow™ Plus
Maximum Aspiration Speed (25mL)	<7 seconds	<5 seconds
Speed Control	Analog	Digital, 6 speeds
Display	n/a	Digital
Dispensing Types	Motor, Gravity and Blowout	
Battery Type	Lithium Ion	
Battery Service Life	8 hours	
Charging Time	2-3 hours	
Pipette Compatibility	1mL to 100mL pipettes, Glass and Plastic	
Internal Filter Type	45µm, hydrophobic	
Weight	200g / 0.44lbs	

ORDERING INFORMATION

Item #	Description	Unit
3385	Motorized serological pipette controller, rechargeable, 120v 60Hz	Each
3385-UK-EU	Motorized serological pipette controller, rechargeable, 240v 50Hz	Each
3387	Motorized serological pipette controller, enhanced, rechargeable, 120v 60Hz	Each
3387-STN	Stand for SeroFlow Plus (3387) serological pipette controller	Each
3387-UK-EU	Motorized serological pipette controller, enhanced, rechargeable, 240v 50Hz	Each
3389-FIL	Replacement filter, 45µm	Each
3389-SILADP	Replacement silicone cone adaptor	Each
3389-BATT	Replacement rechargeable battery	Each
3389--PWRUS	Replacement 120v 60Hz Power Adaptor, with US Plug	Each
3389-PWREU	Replacement 240v 50Hz Power Adaptor, with EU Plug	Each
3389--PWRUK	Replacement 240v 50Hz Power Adaptor, with UK Plug	Each



3385 –
Diamond® SeroFlow™



3387 –
Diamond® SeroFlow™ Plus



www.globescientific.com

Phone 1 (800) 394-4562 • 1 (201) 599-1400
 Fax 1 (201) 599-1406
 E-mail mail@globescientific.com