


12mL Bulb (Whale) Type Urine Centrifuge Tube

Item# 112030

Description

These Globe Scientific "Whale" type urine centrifuge tubes have been designed specifically for urinalysis procedures. The tubes feature a sediment bulb at the bottom and a flared opening which is designed to make pouring off easy.



Features

- 12mL capacity
- 21 x 105mm
- Sediment bulb
- Flared opening
- Polystyrene (PS)
- Graduated to 10mL
- Graduated at 2.5mL, 5mL, 7.5mL and 10mL

Technical Specifications

Spin Rate: 1500 RPM
 Autoclavable: No
 Brittleness temperature: 20°C
 Maximum temperature: 90°C
 Sterility: Non-sterile

Products

Item#	Description	Units
112030	12mL urine centrifuge tube, PS	500 & 1500