

Technical Data Sheet


Rev. 18149

Page 1 of 1

16 x 75mm (8mL Capacity) Test Tubes

Item# 119040





Description

These Globe Scientific test tubes are molded from premium grade polypropylene (PP) for impact strength and stability to common acids, alkalies, and solvents. They can withstand 3000 RPM and are autoclavable and freezable.

Features

- 8mL capacity
- 16 x 75mm
- Polypropylene (PP)
- Non-graduated

Technical Specifications

Spin Rate: 3000 RPM

Autoclavable: Yes

Brittleness temperature: 0°C Maximum temperature: 135°C

Sterility: Non-sterile

Products

Item#	Description	Unit
119040	8mL test tube, with rim, PP	2500



Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406 E-mail mail@globescientific.com Globe Scientific Inc.

400 Corporate Drive • Mahwah NJ 07430

www.globescientific.com