

Snap Lid Multi-Purpose Containers, 8 oz.

Item# 271009



Description

These disposable containers are made of high density polyethylene and are ideal for the collection and storage of histology and pathology specimens, fluids, powders and solids. They are stackable and feature snap on lids. The containers are translucent to aid in viewing of contents.

Features

- Translucent container allows viewing of contents
- Snap-on lid

Technical Specifications


Capacity: 8 oz (240mL)

Sterility: Non-sterile

Material:

Containers: High Density Polyethylene (HDPE)

Lids: Low Density Polyethylene (LDPE)



Products

Item#	Description	Unit
271009	Container, Multi-Purpose, HDPE, 8oz (240mL), Tall Style, Separate Snap Lid, Natural	100