

1mL Nesting Sample Cup, Polystyrene (PS), for 12mm & 13mm Tubes

Item# 5504



Description

This nesting sample cup is used in conjunction with other test tubes or original blood collection tubes when smaller sampling is necessary.

It is ideal for use on direct-tube sampling machines: Olympus® AU480 and AU680, Dade Dimension®, Hitachi® 747, Kodak Ektachem® 700 Series and more.

Features

- 1mL capacity
- For use with 12mm and 13mm tubes
- Excellent clarity

Products

Item#	Description	Unit
5504	1mL nesting sample cup, PS	1000

Technical Specifications

Material: Polystyrene (PS)
 Autoclavable: No
 Temperature range: -10°C to +50°C
 Sterility: Non-sterile

