


.5mL Beckman[®] Type Microcentrifuge Tubes with Attached Snap Cap, Polyethylene (PE)

Item# 110540



Description

These 500uL (0.5mL) Beckman[®] Type microcentrifuge tubes with attached cap are non-graduated. The tubes are made from polyethylene and are not autoclavable. The tubes are sold in bags of 1000.

Features

- Polyethylene (PE)
- Attached snap cap
- Conical bottom
- Non-graduated
- Latex-free

Products

Item#	Description	Packaging
110540	500 uL (.5mL), PE	1000 per case

Technical Specifications

Material: Polyethylene (PE)
 Autoclavable: No
 Temperature range: -50° C to +75°C
 Spin rate: 11,000 RCF
 Shelf life: 5 Years

All instruments and names of manufacturers mentioned herein are the trademarks of their respective companies and are used for reference purposes only. Globe Scientific Inc. expressly disclaims any affiliation to them.