

2.5mL Semi-Micro Spectrophotometer Cuvette, Polymethylmethacrylate (PMMA)

Item# 112147



Description

These Globe Scientific disposable semi-micro cuvettes are designed for use in spectrophotometers. They are produced from UV grade PMMA and include optical windows for optimum transmittance through spectrum range from 300-900nm. These cuvettes are engineered for high reproducibility and offer variations within the range of $\pm 1\%$ absorption. This is vital when batch analysis is being undertaken.

The cuvettes are packaged in dust-free PE foam trays. Each tray contains 100 cuvettes.

Features

- 2.5mL Capacity
- 2 Clear sides
- Excellent clarity

Technical Specifications

Material: UV Grade Polymethylmethacrylate (PMMA)

Dimensions: 12mm x 12mm x 45mm

Path Length: 10mm

Spectrum Range: 300 to 900nm

Autoclavable: No

Sterility: Non-sterile

Products

Item#	Description	Unit
112147	Cuvette, Semi-Micro, 2.5mL, with 2 Clear Sides, UV Grade Polymethylmethacrylate (PMMA)	500