

Graduated Transfer Pipet, 7mL

Item# 135030

Technical Specifications

- Material: Low Density Polyethylene (LDPE)
- Length: 160mm (+/- 0.5mm)
- Capacity: 7mL
- Drops/mL: 20
- Bulb Draw: 3.2mL
- Graduation Increments: 0.5mL
- Graduated to: 3mL
- Sterility: Non-Sterile
- Autoclavable: No
- Temperature Range: -50°C to +75°C

Chemical Resistance

- Excellent resistance (no attack) to dilute and concentrated acids, alcohols, bases and esters
- Good resistance (minor attack) to aldehydes, ketones and vegetable oils
- Limited resistance (moderate attack suitable for short-term use only) to aliphatic and aromatic hydrocarbons, mineral oils, and oxidizing agents
- Poor resistance, and not recommended for use with Halogenated hydrocarbons.

Quality Control

- This product is lot inspected and tested to ensure the highest level of quality.

Recommended Storage Conditions

- Store in cool dry place
- Avoid temperature variations and humidity

