

## Exact Volume Transfer Pipet, 75 $\mu$ L (0.075mL)

Item# 139112

### Technical Specifications

Material: Low Density Polyethylene (LDPE)  
Length: 75mm (+/- 1mm)  
Capacity: 75 $\mu$ L  
Bulb Draw: 75 $\mu$ L  
Sterility: Non-Sterile  
Autoclavable: No  
Temperature Range: -50°C to +75°C

### Chemical Resistance

- Excellent resistance (no attack) to dilute and concentrated acids, alcohols, bases and esters
- Good resistance (minor attack) to aldehydes, ketones and vegetable oils
- Limited resistance (moderate attack suitable for short-term use only) to aliphatic and aromatic hydrocarbons, mineral oils, and oxidizing agents
- Poor resistance, and not recommended for use with halogenated hydrocarbons.

### Quality Control

- This product is lot inspected and tested to ensure the highest level of quality.

### Recommended Storage Conditions

- Store in cool dry place
- Avoid temperature variations and humidity

