


1mL Diamond[®] Essentials[™] Cryogenic Vial, External Thread, Sterile

Item# 3030-1S



Description

Diamond[®] Essentials[™] Cryogenic Vials — a reliable sample storage solution that fits perfectly within your budget. These vials offer exceptional quality without the hefty price tag, making them an ideal choice for researchers seeking top-notch performance.

The range features four volumes that are among the most popular options available in the market today.

Features

- External threaded vial
- Conical bottom, self-standing
- Printed graduations and writing area
- Certified leakproof to 85kPa
- RNase free, DNase free, Human DNA free, Pyrogen free
- Manufactured under the ISO 9001/2015 quality system
- Autoclavable

Packaging:

- 500/case

Products

Item #	Description	Unit
3030-1S	1.0mL, STERILE, external threads, attached screwcap, conical bottom, self-standing, printed graduations, writing space	500

Technical Specifications

Working Volume: 1.0mL

Maximum Volume: 1.0mL

Material: Virgin Polypropylene

Sterility: Irradiation sterilized to SAL10⁻⁶

Temperature Range: -196°C to +121°C

Certified: RNase free, DNase free, Human DNA free, Pyrogen free

Autoclavable: Yes (with cap loosened or removed)