


## 5mL RingSeal™ Cryogenic Vial, External Thread, Round Bottom, Sterile

Item# 3032-5-RB



### Description

RingSeal™ Cryogenic Vials are an excellent choice for long-term sample storage. Built with reliable o-ring seals, this exceptional range offers multiple variants to meet diverse research needs.

Globe Scientific's RingSeal™ Cryogenic Vials are available in both internal and external thread options, as well as self-standing and round bottom formats, and are offered in the most popular volumes.

### Features

- External threaded vial
- Round bottom
- Reliable o-ring seal
- Printed graduations and writing area
- Certified leakproof to 90kPa
- Certified RNase free, DNase free, Pyrogen free
- Autoclavable

### Packaging:

- 500/case

### Products

Item #	Description	Unit
3032-5-RB	5.0mL, STERILE, external threads, attached screwcap, round bottom, printed graduations, writing space	500

### Technical Specifications

Nominal Volume: 5.0mL

Working Volume: 4.5mL

Thread Type: External

Material: Virgin Polypropylene

Sterility: Gamma sterilized to SAL10<sup>-6</sup>

Temperature Range: -196°C to +121°C

Certified: RNase free, DNase free, Pyrogen free

Autoclavable: 121°C for 15 minutes  
(with cap loosened or removed)