

## Cuvette for IL® ACL® Models

Item# 5150



### Description

These Globe Scientific 20-place cuvette rotors are designed for use with the IL® ACL® models: 100, 200, 300, 1000, 2000, 3000, 6000 & 7000.

The cuvette rotors are packaged 25 per sleeve, 4 sleeves per case.

Please note: Not for use on the ACL® 9000 & 10000.

### Features

- 20 Place

### Technical Specifications

Material: Polystyrene (PS)

Autoclavable: No

Sterility: Non-sterile

### Products

Item#	Description	Unit
5150	Rotor Cuvette, 20-Place, Polystyrene (PS) 25 per sleeve, 4 sleeves per case	100

*Note: All instruments and names of manufacturers mentioned herein are the trademarks of their respective companies and are used for reference purposes only. Globe Scientific Inc. expressly disclaims any affiliation to them.*