

12mL Conical Bottom Centrifuge Tube

Item# 112040

Description

These Globe Scientific general purpose centrifuge tubes feature conical bottoms and are made from polypropylene (PP). The tubes are graduated at 0.5, 1, 2, 5 and 10mL.

Features

- 12mL capacity
- 16 x 100mm
- Conical bottom
- Polypropylene (PP)
- Graduated to 10mL
- Graduated at 0.5mL, 1mL, 2mL, 5mL and 10mL

Technical Specifications

Spin Rate: 3000 RCF

Sterility: Non-sterile

Products

| Item# | Description | Units |
|--------|---|------------|
| 112040 | 12mL centrifuge tube, PP, 500/bag, 4 bags/case | 500 & 2000 |

