


## 1-300µL Filter Pipette Tips, Low Retention, STERILE

Item# 150925



### Technical Specifications

Material: Polypropylene (PP)  
 Autoclavable: Yes  
 Maximum Volume: 300µL  
 Sterility: Gamma Irradiation SAL10<sup>-6</sup>

### Products

Item #	Description	Unit
150925	Filter Pipette Tip, 1-300uL, Certified, Universal, Low Retention, 59mm, Natural, STERILE, 96/Rack, 10 Racks/Box, 2 Boxes/Carton	1920

### Description

These Globe Scientific STERILE filter pipette tips are manufactured from clear high-quality polypropylene (PP). Special processing renders these tips hydrophobic for exceptionally low sample retention. The tips are certified DNase free, RNase free and non-pyrogenic. These premium quality pipette tips are manufactured in a Class 100,000 cleanroom environment under the ISO13485:2016 and ISO9001:2015 quality management systems. A universal fit design provides an airtight fit on most brands of pipettors. The tips are offered in extra-rigid autoclavable racks.

### Features

- Manufactured from clear polypropylene (PP)
- Polypropylene (PP) filters protect samples and pipettors from contamination
- DNase-free, RNase-free and non-pyrogenic
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485:2016 and ISO9001:2015 quality management system
- Sterilized by gamma irradiation SAL10<sup>-6</sup> (ISO11137)
- Packaged in extra-rigid autoclavable racks
- Racks are printed with lot number for quality traceability

### Material:

**Tips:** Polypropylene (PP), Color: Natural

**Filters:** Polypropylene (PP), Color: White

**Rack Box:** Polypropylene (PP), Color: Grey/Red

### Packaging:

- 96/Rack, 10 Racks/Box, 2 Boxes/Carton

### Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400  
 Fax 1 (201) 599-1406  
 E-mail mail@globescientific.com

### Globe Scientific Inc.

400 Corporate Drive • Mahwah, NJ 07430  
 www.globescientific.com