

## Product **Spec Sheet**

**Globe Item#:** V-3015T

**Description:**

Cryogenic Vial, 5.0mL, STERILE, External Threads, Attached Screw Cap, Self-Standing, No Printed Graduations or Marking Place, Bulk

**Material:** Polypropylene (PP), USP Class VI


**Specifications:**

Gamma Sterile  
RNase/DNase Free  
Cytotoxin Free, Non-Pyrogenic  
IATA Compliant (95kPA)

**Dimensions:**

12mm x 85mm without Cap  
12mm x 91mm with Cap

**Packaging:** Case of 500, Bulk, Sterile



**Customer service, sales and technical support:**

Phone 1 (800) 394-4562 • 1 (201) 599-1400  
Fax 1 (201) 599-1406  
E-mail [mail@globescientific.com](mailto:mail@globescientific.com)

**Globe Scientific Inc.**

400 Corporate Drive • Mahwah, NJ 07430  
[www.globescientific.com](http://www.globescientific.com)