

Weighing Boats & Dishes


A selection of premium weighing boats and dishes for a wide range of applications.



Weighing Boats & Dishes

Globe Scientific Weighing Boats are made for a wide range of applications including weighing, dispensing, and storing liquids, powders, and solids. Weighing boats can also be used as quick freeze trays for sample material, discard trays for broken ampules, or mixing trays for small batches. The weighing boats are lightweight and durable with smooth surfaces and sloping sides to facilitate accurate pour-outs with no sample loss, and flat bottoms that prevent tipping. Made from polystyrene, the weighing boats are engineered to flex for controlled transfers. The weighing boats are available in hexagonal, square, diamond and pour boat styles, all in a variety of sizes.

- Made from polystyrene
- Engineered to flex for controlled pouring
- Biologically inert, contaminant-free, and resistant to most dilute acids, alcohols, bases, and aqueous solutions
- Food-grade compliant
- Withstand temperatures -10C to 80°C
- Hydrophobic surface prevents absorption of sample
- Anti-static, non-stick surface prevents sample loss
- Flat bottoms ensure stability and prevent tipping
- Smooth surfaces and sloping sides facilitate accurate sample transfers
- Stackable for easy storing
- Disposable for quick and easy clean up



Hexagonal Weighing Boats

- Flexible polystyrene enables the weigh boats to bend easily to form pour spouts for controlled transfers

Ordering Information



ESD safe

Item #	Description	Unit
3615	20mL, Anti-static, Flat Bottom	500
3616	50mL, Anti-static, Flat Bottom	500
3617	200mL, Anti-static, Flat Bottom	500
3618	350mL, Anti-static, Flat Bottom	500



Diamond Shaped Weighing Boats

- Flexible polystyrene enables the weigh boats to bend easily to form pour spouts for controlled transfers

Ordering Information



ESD safe

Item #	Description	Unit
3614-5	5mL, 58 x 36 x 6mm, Anti-static	250
3614-30	30mL, 80 x 60 x 18mm, Anti-static	250
3614-100	100mL, 125 x 100 x 18mm, Anti-static	250

Weighing boats are also ideal for applications outside the lab, such as measuring and weighing food, digital scale trays, weighing and mixing dyes and epoxy, measuring powders and liquids, disposable dishes for glues and paints, weighing medications, and more.





Pour Boat Weighing Dishes

- Circular OR square flat bottom options ensure stability and prevent tipping

Ordering Information



ESD safe

Item #	Description	Unit
3623	10mL, 57 x 44 x 8mm, Anti-static, Circular Flat Bottom	250
3624	140mL, 137 x 89 x 25mm, Anti-static, Circular Flat Bottom	250
3625	270mL, 191 x 121 x 25mm, Anti-static, Circular Flat Bottom	250
3626	15ml, 53 x 38 x 13mm, Anti-static, Square Flat Bottom	250
3627	25ml, 68 x 50 x 16mm, Anti-static, Square Flat Bottom	250
3628	120ml, 133 x 85 x 26mm, Anti-static, Square Flat Bottom	250
3629	240ml, 181 x 110 x 26mm, Anti-static, Square Flat Bottom	250

Square Weighing Boats

- Circular OR square flat bottom options ensure stability and prevent tipping
- White and black color options available to contrast with samples colors for improved sample recovery

Ordering Information



ESD safe

Item #	Description	Unit
3620	7mL, 41 x 41 x 8mm, Anti-static, Circular Flat Bottom, White	500
3621	100mL, 89 x 89 x 25mm, Anti-static, Circular Flat Bottom, White	500
3622	330mL, 140 x 140 x 25mm, Anti-static, Circular Flat Bottom, White	500
3619-10	10ml, 45 x 45 x 8mm, Anti-static, Square Flat Bottom, White	500
3619-100	100ml, 80 x 80 x 24mm, Anti-static, Square Flat Bottom, White	500
3619-250	250ml, 135 x 135 x 24mm, Anti-static, Square Flat Bottom, White	500
3619-10BK	10ml, 45 x 45 x 8mm, Anti-static, Square Flat Bottom, Black	500
3619-100BK	100ml, 80 x 80 x 24mm, Anti-static, Square Flat Bottom, Black	500
3619-250BK	250ml, 135 x 135 x 24mm, Anti-static, Square Flat Bottom, Black	500



Pour boat weighing dishes available in 2 styles.



Square weighing boats with circular bottoms.



Square weighing boats with square bottoms offered in white and black.



Aluminum Weighing Dishes

A variety of styles that can be used for applications including general purpose weighing, dispensing, storage, balance pans, food processing, moisture testing, embedding mold, evaporating, dust covers and more.

- Made from lightweight sturdy aluminum
- Disposable
- Oil residue free and non-coated
- Smooth flat bottoms ensure stability and prevent tipping
- Stackable for easy storing and angled sides enable easy separation when ready for use



Aluminum weighing dishes with crimped walls and finger grip tab.



Aluminum weighing dishes with crimped walls and curled lip.

Aluminum weighing dishes with smooth walls. 100mm dish has flanged lip.

Aluminum weighing/drying pan with flanged lip.

Ordering Information

Item #	Description	Unit
8305	20mL, 50mm OD, with Tab Crimped Walls, No Lip	100
8306	40mL, 60mm OD, with Tab Crimped Walls, No Lip	100
8307	80mL, 76mm OD, with Tab Crimped Walls, No Lip	100
8309	8mL, 30mm OD, with Tab Crimped Walls, No Lip	100
8312	20mL, 51mm OD, Crimped Walls, with Lip	100
8313	75mL, 70mm OD, Crimped Walls, with Lip	100
8314	70mL, 70mm OD, Crimped Walls, with Lip	100
8315	200mL, 127mm OD, Crimped Walls, with Lip	50
8308	70mL, 72mm OD, Smooth Walls, No Lip	100
8310	150mL, 100mm OD, Smooth Walls, Flanged Lip	50
8323	Pan, 60mL, 102mm OD, Smooth Walls, Flanged Lip	80

Micro Aluminum Weighing Dishes

- Made from lightweight sturdy aluminum
- Disposable
- Oil residue free and non-coated
- Smooth flat bottoms ensure stability and prevent tipping
- Stackable for easy storing and angled sides enable easy separation when ready for use
- Packaged in easy-to-use plastic containers for storage and dispensing
- Crimped walls for rigidity
- Smaller surface area reduces evaporation rate
- Choose between round or oval weighing dishes



Ideal for microbalances, organic micro elementary and thermal analysis, and pharmaceutical use.

Ordering Information

Item #	Description	Unit
8326	1.12mL Micro, 25L x 7W x 7H mm, Rectangle, Crimped Walls, No Lip	100
8327	1.12mL Micro, 25L x 7W x 7H mm, Rectangle, with Tab, Crimped Walls, No Lip	100
8328	.35mL Micro, 13mm OD x 3.5H mm, Round, with Tab, Crimped Walls, No Lip	100
8329	10mL Micro, 65L x 18W x 9H mm, U Round, Crimped Walls, No Lip	100



www.globescientific.com

Globe Scientific Inc.

400 Corporate Drive • Mahwah, New Jersey 07430 • USA

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com